
# User Manual

Please read the manual before using the product.

## Contents

Safety Precautions .....	4
Product Introduction .....	4
Product Specifications .....	5
Product Description .....	5
How To Use .....	6-8
Troubleshooting .....	8
安全上のご注意 .....	9
製品説明 .....	9
製品仕様 .....	10
製品構成 .....	10-11
使用方法 .....	11-13
ユーザーサポート及びQ&A .....	13
Warranty .....	14

## Safety Precautions

### Warning

1. If the device is disassembled by the user, no replacement or refund is provided.
2. Do not stare directly into the laser. This may cause blindness.
3. The device should not be used by pregnant women and children under the age of 13.
4. The recommended usage time is 30 minutes per session, up to 3 sessions per day. May cause burns, skin irritation and allergic reactions to sensitive skin.
5. Patients with cardiovascular disease, kidney disease, cancer, or diabetes should not use this device.
6. Keloid patients and anti-inflammatory drug users should not use this device.
7. Patients with hypersensitive skins, hematuria, or hypertrophic scar marks should not use this device.
8. Do not use if you have sensitive skin with allergies. It may cause skin irritation.
9. Do not use with other medical devices or electronics such as wearable electronic medical devices, life support electronics, heartbeat aids, electromagnetic generating devices, etc.
10. This device contains a built-in battery. External pressure, impact, heat or fire may cause explosion. For example, do not place this device close to or inside of a stove, microwave, heat cooker or high-pressure container.

### Caution


1. Please use and store in a dry place, without water and moisture.
2. Avoid external vibration or impacts when storing, using or carrying the device.
3. Read the Usual Manual in case of damage, malfunction, or non-operation.
4. Attach the laser surface tightly onto the body area for complete contact with the device.
5. Do not use on deep wounds or bleeding wounds.
6. This is multiple wavelength infrared light therapy device manufactured by WellsCare. Use the product only for the purposes that it is intended for. Use the device with the genuine
7. Do not force the device into the charger's port for recharge. It may damage the charging port.
8. Using charging cables other than the one included in the package may cause damages to the product.
9. Do not use the device with cosmetics, cream, gel, or ointment that are not manufactured by WellsCare.
10. The heat you may feel from the device is not produced by the battery but is the normal heating function of the device.

### Liability

- We do not take any responsibility for any problem or damage caused by user's carelessness.
- We do not take any responsibility for any problems or damage caused by the user disassembling and modifying the device.

## Product Introduction

Low Level Laser, also called LLL or cold laser, refers to laser with 500mW wavelength or lower. Based on various clinical research and studies, the low level laser therapy (LLLT) is known to help with both chronic and acute pain relief and treatment.

Especially, the 600nm-1000nm laser wavelength effectively works for pain relief, blood flow improvement, vasodilation, lymph node stimulation, cell activation by its biochemistry mechanism.

Our product is a compact wearable laser device that emits dual wavelengths, 650nm and 830nm. It is developed as an everyday home care cold laser device that anyone can use as per their convenience to shoulders. This mechanism is known to promotes generating ATP(Adenosine Triphosphate) cells that make energy from the body and protein.


※ This is a healthcare massage device with investigate wavelength infrared irradiation and vibration functions.

## Product Specifications

Product Name	Dual Wavelength Infrared Light Device
Model Name	IASO-SP1
Manufacturer	WellsCare.Co.Ltd
Exterior Size	48mm x 48mm x 12mm
Packaging Unit	1 Set
Laser Medium	650nm(InGaAlP) / 830nm(GaAs)
Laser Class	Class 3R (Low-level laser)
Laser Power	Less than 5.0mW
Output Type	Continuous
Emmission Time	30min ± 10%
Vibrator	Less than 15,000 rpm (Fixed type)
Battery	3.7V, 260mAh, 0.962Wh (Li-Polymer)
Weight	28g
Charging Port	5V/2A micro USB Type
Rated input	input power 5.0V, 2,000mA (micro USB Type Cable)
Charging Duration	<p><b>Charging time:</b> approximately 60 minutes for a full charge</p> <p><b>Operation Time :</b> [When Mode 1 is in use] 2 x 30-minute session (Approximately 60 minutes)                      [When Mode 2 or 3 is in use] 30-minute session + 20-minute session (Approximately 50 minutes)</p>
Storage method	Avoid hot, damp, unclean, and dusty areas when storing the device.

## Product Description

### 1. Main Device

No	Category	Description
1	Main Body	Parts for operation of the device and laser movement are present inside the device.
2	Switch, Operation Display	Located on the side of the instrument, controls the on/off switch of the device. Operation display: Red for charging, blue for product operation, green for laser output. * Blue light turns on when fully charged.
3	Charging Port	A port to recharge the main body of the device.
4	Laser Output Port	A port for laser emission Center laser output port emits 830 nm wavelength laser (invisible) Surrounding 4 output ports emit 650 nm wavelength laser.


#### Caution

- For the first time use after purchase, charge for a sufficient time before use.
- Surrounding 4 output ports emit 650 nm wavelength laser.
- Warranty is not provided for damages caused by incorrect use of the charger.
- Charge for over an hour, if the device has not been use for more than a week.

## 2. Charging Cradle

Image	 charging port	 charging port
Name	<b>2 Channel Charging Cradle</b> (For Double Package Only)	<b>4 Channel Charging Cradle</b> (For Quadruple Package only)
Purpose	Simultaneously Charges up to 2 devices	Simultaneously Charges up to 4 devices
Size	72 x 84 x 33 (mm)	72 x 104 x 33 (mm)
Rating	DC 5V,2A (micro USB 5pin)	

## 3. Accessories

Hydrocolloid  
Adhesive Patches


Band

Portable Tray


Charging Cable


User Manual

## How To Use

### 1. Main Device


Press and Holding the button on the side for 2 seconds powers up the device with an alarm sound. A blue light will be displayed to indicate standby mode. Press and Holding the button again for 2 seconds turns off the device.

The device starts operating and its modes can be changed by pressing the button on the side while the device is on (blue light).

Caution

The button may not operate if it is pressed from the top.


#### Mode 1 (Laser Only)

Pressing the button once from the Stand-by mode, will initiate the laser only function.

#### Mode 2 (Laser + Continuous Vibration)

Pressing the button twice from the Stand-by mode, will initiate the laser with continuous vibration function.

#### Mode 3 (Laser + Pulse Vibration)


Pressing the button Three times from the Stand-by mode, will initiate the laser with pulse vibration function.

## When the battery is low



- When the battery is less than 10%, the device vibrates with flashing LED light. This is not an operation error. Simply turn off the device and recharge it.
- When charging, red LED light turns on. Start using the device when it turns to green, which means that the device is fully charged.

## 2. How to use the portable tray


Place the device inside the portable tray with the button aligned to the lower part of the portable tray. The button may get damaged internally if the device is placed with the button facing against the higher part of the portable tray.

[Correct use]

[Incorrect use]


## 3. How to use the activity band

Insert the device from the back of the activity band with the button aligned to the narrower part of the device insert on the band. Slide the device first into the narrower part of the band and fasten.



The button may get damaged internally if it is not correctly placed into the narrower part of the device insert on the band.

[Correct use]

[Incorrect use]


### How to insert the device after wearing the band

Make sure to slide the charging port area of the device into the band first after placing the band around your knee or wrist. Then, carefully fasten the button area of the device into the other side of the band.


It is recommended to wear the activity band first after adjusting its length to fit your knee or wrist before inserting the device.


[On Knees]

[On Wrists]


## 4. How to charge the device



### Correct use of 2-port charging cradle

Insert the devices in the charging cradle with the buttons facing upward. When using a 2-port charging cradle, make sure that the interior sides of the devices are facing each other.

### Correct use of 4-port charging cradle

Insert the devices in the charging cradle with the buttons facing upward. When using a 4-port charging cradle, make sure to place all devices facing away from the charging cable.

### When charging the device directly with a charging cable


Make sure to check the shape of the cable and the charging cradle port before inserting the cable into the port.

### Why won't the device charge?


- Make sure that the pins inside the charging cradle are intact.
- Make sure to insert the 5-pin charging cable into the charging port of the cradle to the correct direction.
- Make sure to insert the devices in the charging cradle in the correct position.

## 5. How to use Adhesive Patches


The patches are made of hydrocolloid material and do not cause irritation for most skin types. However, if you are allergic to hydrocolloid, refrain from using the patches.

Open the pouch and take out a sheet of adhesive patches.

Attach the adhesive patch, carefully aligning it with the device starting from the red part of the patch.

Remove the protective film by peeling it off from the small protruding edge. Then, attach the device on the desired place on your body.

### Note while using the Adhesive Patches

- Please note that adhesivity of the patches may deteriorate if exposed for a prolonged time.
- Make sure to clean and dry the skin surface before attaching the device for the best result.
- Do not attach a patch on top of a used patch. Remove the used patch from the device before attaching a new one.

## Troubleshooting

Please read the manual thoroughly before use. Please follow all the steps above before contacting the service center.

Symptoms	Solutions
The device does not turn on.	<ol style="list-style-type: none"> <li>1. Make sure that the battery is charged.</li> <li>2. Press and hold the power button for more than 2 seconds.</li> </ol>
The device is on but does not operate.	<ol style="list-style-type: none"> <li>1. The device does not operate while charging.</li> <li>2. Make sure to turn on the device by pressing the button and the green light is on. * A blue light indicates the Stand-by mode.</li> </ol>
The device does not turn on even after charging for a long time.	<ol style="list-style-type: none"> <li>1. Make sure that the connector of the charger is properly connected to the device and charge again.</li> <li>2. Clean the charger port and make sure that there is no foreign objects in the port.</li> <li>3. The charger may be defective. Try charging the device with a different charger. * A strong impact on the charging cradle may damage the charging pins inside the cradle.</li> </ol>
The device gets hot while in use.	<ol style="list-style-type: none"> <li>1. The device may produce heat while in use and it is normal.</li> <li>2. When the heat feels too hot or causes discomfort, reduce the usage time per session.</li> </ol>
The device is not emitting laser light.	<ol style="list-style-type: none"> <li>1. The laser emitting from the center of the device is almost infrared and is invisible.</li> <li>2. Due to the characteristics of the product, the laser sensitivity and brightness may vary and it is normal.</li> </ol>
Vibration is not strong enough or does not work.	<ol style="list-style-type: none"> <li>1. Vibration intensity may feel different depending on the user. Different users may feel more or less intense vibration.</li> <li>2. When vibration does not work, restart the device by turning off the device and on again.</li> </ol>
The device is making a noise.	<ol style="list-style-type: none"> <li>1. It is normal that the device produces minor noises while vibrating.</li> <li>2. The vibration noises may sound louder or less loud depending on the environment.</li> </ol>

### Exchange and Return

- Replacement or return available within 14 days after purchase.
- No replacement or return after the packaging box is opened.
- No replacement or return for damage caused by the user or loss of the product.
- Contact the service center if the device is still not working after checking the symptoms and solutions listed above.

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
**Date of purchase :**

### **Detail**

1. Warranty period : 1 year
2. Paid maintenance/repair is available after the warranty expires.
3. Product parts will be available for up to a year after product discontinuation.

### **Exclusions**

- Break down of the product caused by user error/injury
- In case of damage/failure caused by external shocks
- If the product is used for purposes other than its original purpose
- If the user arbitrarily disassembled and modified the product
- Corrosion due to water ingress of foreign substances contamination
- Parts/supplies/accessories etc. Not supplied by the manufacturer
- Other user's fault is recognized

**Manufactured by :**

WellsCare Co., Ltd.

507, 107, Gwanggyo-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do,  
Republic of Korea

**Distributed by**

WellsCare Inc.

3580 WILSHIRE BLVD, STE 1460, LOS ANGELES, CA 90010, US

Website : [www.wellscare.net](http://www.wellscare.net) / [www.wellscare.kr](http://www.wellscare.kr)

E-mail : [support@wellscare.net](mailto:support@wellscare.net)

Tel : (213) 205-7051

[www.wellscare.net](http://www.wellscare.net)